

# SURGICAL MASKS

- BREATHABLE
- COMFORT
- FIT DESIGN
- SPLASH RESISTANT

50 pcs



**EAR LOOP  
(3 ply)**



**HEAD LOOP  
(3 ply)**



**TIE ON  
(3 ply)**



## Made From High Quality Materials :

- Melt-blown nonwoven (MB.25 gsm\*175 mm)
- White Non-wovens (26 gsm\*200 mm)
- Elastic Ribbon (33 mm) / Earloop Tie
- Blue Non-wovens (28 gsm\*175 mm)
- Plastic Coated Strip