

SH 210

Portable Ventilator

ETERNITY



CE 0123

Technical specification

Physical Specification

Display: 5.6" TFT color Screen

Dimension: 480mm (H) * 260mm (W) * 420mm (D)

Weight: 7kg (Gross Weight)

Ventilation Specification

Ventilation Mode : A/C ,SIMV, SPONT, SIGH, MAN

VT : 50~1500 mL

Frequency 4~80 bpm

I: E 3:1~1:5

Inspiratory time 0.2~10s

Pressure trigger sensitivity

Flow trigger sensitivity 1~20 L/min

PEEP 0~20 cmH₂O

Pinsp 5~80 cmH₂O

FiO₂ 45~100%

Mechanical PEEP (Option)

Monitoring

VTI, VTE, MV, FIO₂, Frequency, Fspn,

PEEP, Ppeak, Pmean

Waveform : P-T, F-T, V-T,

Loops : P-V loop, V-F loop

Alarm

Frequency Upper & Low

Continuous airway pressure high

Apena, O₂ failure

AC Power Source failure, Low battery

Battery exhausted

Alarm silence time: ≤120 S

VTE Low/ High

MV Low/ High

FIO₂ Low/ High

Paw Low/ High

Optional Configuration

Trolley Breathing Circuit Arm

Carry Bag Humidifier

Carry Steel

BEIJING ETERNITY ELECTRONIC TECHNOLOGY ., LTD

Add : No .17 Xijing Road , Badachu , Zhonguancun High-tech Zone

Shijingshan District , Beijing 100041 P.R China

TEL: 0086 - 10 - 68827491

FAX : 0086 - 10 - 68817023

www.ernitymed.com

The logo for ETERNITY features the word "ETERNITY" in a bold, blue, sans-serif font. A red, stylized crescent shape is positioned above the letter "N", resembling a rising sun or a medical symbol.