

Robust. Powerful. Efficient.

KaVo **expert** series

KaVo EXPERTtorque LUX E680 L

- **Max. working efficiency**, high torque and optimal power of up to 18 W
- **Optimal service life** due to its excellent head-protect system and stainless steel heads
- **Proven configuration efficiency** due to customer-replaceable water filter and robust chuck system

KaVo EXPERTtorque Mini LUX E677 L

- **Perfect vision** due to smaller head size and 100°/19° head/elbow angle combination
- **Efficiency** due to robust KaVo chuck system for standard and short-stem drills
- **Optimal service life** due to high chuck retention force of up to 30 N



KaVo EXPERTtorque
LUX E680 L

KaVo EXPERTtorque
Mini LUX E677 L



KaVo EXPERTtorque Turbines

KaVo EXPERTtorque LUX E680 L

Power:

- 18 watts power
- Stable speed

Improved vision:

- Head height reduced to 13.6 mm
- 100°/19° head/elbow angle combination (KaVo patent)
- Integrated light

Efficiency and infection protection:

- Four-nozzle spray
- Head-protect system: New improved back-suction effect while the turbine comes to a stop
- Chuck system: Low trigger force

Low noise level:

- Up to 62 dB(A)

Better hygiene conditions:

- Spray air/spray water guided on the interior
- New customer-replaceable water filter



KaVo EXPERTtorque Mini LUX E677 L

Improved vision:

- 12.2 mm head height (miniature head)
- 100°/19° head/elbow angle combination (KaVo patent)
- Integrated light

Power:

- 15 watts power
- Stable speed

Efficiency and infection protection:

- Three-nozzle spray
- Head-protect system: New improved back-suction effect while the turbine comes to a stop
- Chuck system: Low trigger force, use of standard and short-stem drills

Low noise level:

- Up to 65 dB(A)

Better hygiene conditions:

- Spray air/spray water guided on the interior
- New customer-replaceable water filter



Special features of the EXPERTtorque:

- EXPERTgrip surface for a safe grip
- Customer-replaceable water filter for consistently good spray
- Chuck retention force of up to 30 N
- Carbide guide sleeve for perfectly running drill
- Original KaVo ceramic ball bearing
- 18-month warranty
-  Made in Germany

Handpiece servicing:

KaVo QUATTROcare PLUS 2124 A

- Optimal servicing quality with QUATTROcare Plus Spray
- Simple and comfortable handling and rapid programme flow
- Low initial purchase and consumable costs



KaVo expert series

| Article | Description | Material number |
|-------------------------|---|-----------------|
| E680 L | EXPERTtorque LUX | 1.006.8700 |
| E680 C | EXPERTtorque (without light) | 1.006.9000 |
| E680 LS | EXPERTtorque LUX, for Sirona couplings | 1.006.4600 |
| E680 LN | EXPERTtorque LUX, for NSK couplings | 1.006.4700 |
| E677 L | EXPERTtorque Mini LUX | 1.007.3600 |
| E677 LN | EXPERTtorque Mini LUX, for NSK couplings | 1.006.0100 |
| MULTIflex 465 LED | MULTIflex LUX coupling 465 LED | 1.007.3197 |
| MULTIflex 465 RN | MULTIflex coupling 465 RN – without light | 0.553.1600 |
| QUATTROcare PLUS 2124 A | Handpiece servicing system | 1.008.3805 |



KaVo. Dental Excellence.