



# 1. Main body

- Maximum load up to 185kgs
- The color of the unit chair is not binding according to the production at the time or special request
- Integrated Trendelenburg movement for patient comfort
- Extended positions:
  - o Highest - 830mm
  - o Lowest - 350mm



1.

# 2. Your practice in the best light

- Natural white light for every selection of color shades
- COMPOsave mode for delayed hardening of the filling material



2.

# 3. It's all about your back

- Ergonomics and comfort in perfect harmony
- Healthy sitting posture have a significant influence on the profitability of your practice



3.

# 4. Prevention - Hygiene

- Daily hygiene more efficient with easy to clean surface
- Automated assistant Hygiene unit
- Protect water lines from contaminated water with KaVo OXYGENAL 6



4.



## 5. Show it to your patients

- KaVo Screen HD – 19 inch medical grade monitor, for hygiene safety:
  - o Safety glass panels
  - o Narrow gap on the housing for easy cleaning

5.



## 6. Efficiency for your assistant

- Sensitive touch panel for quick and easy function control
- Flexible height adjustment and extensive swivel range, pure ergonomics


6.



## 7. Type Table

- TM Table
- The basics within your reach
- Ergonomics – makes your daily work simple

7.



## 8. Connects what belongs together

- More efficiency in your everyday routine
- Screenshot option of a region of interest and show the patient
- Saves time and no need additional processing
- Every diagnosis is logged, summarized and saved completely

8.

"Every treatment is different,  
but in every moment I can  
rely on perfect ergonomic."

Dr. Piotr Wisniewski  
Dentist, Głuchołazy

## Accessories

### KaVo EXPERTtorque LUX E680 L

- **Max. working efficiency**, high torque and optimal power of up to 18 W
- **Optimal service life** due to its excellent head-protect system and stainless steel heads
- **Proven configuration efficiency** due to customer-replaceable water filter and robust chuck system



### KaVo EXPERTtorque Mini LUX E677 L

- **Perfect vision** due to smaller head size and 100°/19° head/elbow angle combination
- **Efficiency** due to robust KaVo chuck system for standard and short-stem drills
- **Optimal service life** due to high chuck retention force of up to 30 N



### KaVo EXPERTmatic series

#### KaVo EXPERTmatic LUX E20 L

- Transmission ratio 1:1
- Easy access due to a small head
- Efficient design with single spray and internal spray media
- E20 L with glass rod light conductor,



#### KaVo EXPERTmatic E10 C

- Handpiece, transmission ratio 1:1
- Efficient design with single spray and internal spray media
- Suitable for handpiece and contra-angle drill

