

KaVo Primus™ 1058 Life TM

Lighting

KaVoLux 540 LED Operating Light

Dentist Element

1x Turbine with 465 LED
Coupling
1x703 LED
1x PiezoLED Scaler
1x Triple Function Syringe
1x ERGOcam One

Assistant Element

1x Saliva Injector
1x Hi Vac Suction
1x Triple Function Syringe

Accessories

1x Standard Foot Control
1x 19 Inch Medical Grade
HD Monitor
2x KaVo PHYSIO Evo Stool
1x Patient Service Box

Hygiene Function

1x DVGW Water Kit

Handpieces

4x EXPERT Series
E680L
E677L
E20L
E10C

Type Table

Compact Chair : TM Table



Made in Germany

KaVo. My Family of Choice

ESTETICA
E30

PRIMUS
1058 LIFE

ESTETICA
E50 LIFE

ESTETICA
E70/E80 VISION

KAVO
Dental Excellence

1. Main body

- Maximum load up to 185kgs
- The color of the unit chair is not binding according to the production at the time or special request
- Integrated Trendelenburg movement for patient comfort
- Extended positions:
 - o Highest - 830mm
 - o Lowest - 350mm

1.



2. Your practice in the best light

- Natural white light for every selection of color shades
- COMPOsave mode for delayed hardening of the filling material

2.



3. It's all about your back

- Ergonomics and comfort in perfect harmony
- Healthy sitting posture have a significant influence on the profitability of your practice

3.



4. Prevention - Hygiene

- Daily hygiene more efficient with easy to clean surface
- Automated assistant Hygiene unit
- Protect water lines from contaminated water with KaVo OXYGENAL 6

4.





5. Show it to your patients

- KaVo Screen HD – 19 inch medical grade monitor, for hygiene safety:
 - o Safety glass panels
 - o Narrow gap on the housing for easy cleaning

5.



6. Efficiency for your assistant

- Sensitive touch panel for quick and easy function control
- Flexible height adjustment and extensive swivel range, pure ergonomics


6.



7. Type Table

- TM Table
- The basics within your reach
- Ergonomics – makes your daily work simple

7.



8. Connects what belongs together

- More efficiency in your everyday routine
- Screenshot option of a region of interest and show the patient
- Saves time and no need additional processing
- Every diagnosis is logged, summarized and saved completely

8.

"Every treatment is different,
but in every moment I can
rely on perfect ergonomic."

Dr. Piotr Wisniewski
Dentist, Głucholazy

Accessories

KaVo EXPERTtorque LUX E680 L

- **Max. working efficiency**, high torque and optimal power of up to 18 W
- **Optimal service life** due to its excellent head-protect system and stainless steel heads
- **Proven configuration efficiency** due to customer-replaceable water filter and robust chuck system



KaVo EXPERTtorque Mini LUX E677 L

- **Perfect vision** due to smaller head size and 100°/19° head/elbow angle combination
- **Efficiency** due to robust KaVo chuck system for standard and short-stem drills
- **Optimal service life** due to high chuck retention force of up to 30 N



KaVo EXPERTmatic series

KaVo EXPERTmatic LUX E20 L

- Transmission ratio 1:1
- Easy access due to a small head
- Efficient design with single spray and internal spray media
- E20 L with glass rod light conductor,



KaVo EXPERTmatic E10 C

- Handpiece, transmission ratio 1:1
- Efficient design with single spray and internal spray media
- Suitable for handpiece and contra-angle drill

