



CaviCide™

Surface Disinfectant/Decontaminant Cleaner

A convenient, ready-to-use, intermediate-level disinfectant.

- Effective against TB in 3 minutes, and HIV-1, HCV, HBV and MRSA in 2 minutes*
- Disinfectant, cleaner and decontaminant all in one
- For use on hard, nonporous surfaces, even in neonatal and infant care areas

Part No.	Product Description	Packaging
13-1000	CaviCide – 1 gallon	4 bottles/case

Composition/Information on Ingredients

Component	CAS No.	Amount
Water	7732-18-5	70-80%
Isopropanol	67-63-0	17.2%
Ethylene Glycol Monobutyl Ether (2-Butoxyethanol)	111-76-2	1-5%
Diisobutylphenoxyethoxyethyl dimethylbenzyl ammonium chloride	121-54-0	0.28%

Handling and Storage

Precautions for Safe Handling: Avoid contact with eyes, skin and clothing. Wear appropriate eye protection when handling. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse. Combustible liquid. Keep away from heat, open flames and all other sources of ignition. Do not smoke in storage or use areas. Keep containers closed when not in use. Do not reuse empty containers.

Empty containers retain product residues and may be hazardous. Do not flame cut, drill, weld, etc. on or near empty containers, even empty.

Conditions for Safe Storage, Including any Incompatibilities: Store in a cool, well ventilated area away from heat, oxidizers and all sources of ignition. Do not contaminate water, food or feed by storage.