



CaviCide1™

1-Minute, 1-Step Surface Disinfectant Liquid

Fast acting, intermediate-level surface disinfectant spray available in many sizes to meet a variety of needs, including cost conscious and large volume users.

- Effective against TB, Norovirus, *A. baumannii* and other microorganisms in just 1 minute*
- 1-Step cleaning and disinfecting allows for faster room turnover
- Low alcohol formula offers excellent materials compatibility
- For use on hard, nonporous surfaces, even in neonatal and infant care areas

Part No.	Product Description	Packaging
13-5000	CaviCide1 - 1 gallon	4 bottles/case

Composition/Information on Ingredients

Component	CAS No.	Amount
Water	7732-18-5	70-80%
Isopropanol	67-63-0	15%
Ethanol	64-17-5	7.5%
Ethylene Glycol Monobutyl Ether (2-Butoxyethanol)	111-76-2	1-5%
Didecyldimethylammonium chloride	7173-51-5	0.76%

Handling and Storage

Precautions for Safe Handling: Avoid contact with the eyes. Avoid prolonged contact with skin and clothing. Use with adequate ventilation. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Flammable liquid and vapor. Keep away from heat, sparks, open flames and all other sources of ignition. Do not smoke in storage or use areas. Keep containers closed when not in use. Do not reuse empty containers.

Empty containers retain product residues and may be hazardous. Do not flame cut, drill, weld, etc. on or near empty containers, even empty.

Conditions for Safe Storage, Including any Incompatibilities: Store in a cool, well ventilated area away from heat, oxidizers and all sources of ignition.