

CaviWipes™

Disinfecting Towelettes

A convenient, ready-to-use, intermediate-level disinfecting wipe available in three sizes.

- Kills TB in 3 minutes and MRSA, HIV-1, HBV and HCV in 2 minutes*
- CaviCide and Towelette in one
- For use on hard, nonporous surfaces
- Sturdy and durable, nonwoven towelette will not bunch up during use



Part No.	Product Description	Packaging
13-1100	CaviWipes (160 wipes per can)	12 canisters/case
13-1150	CaviWipesXL (65 wipes per can)	12 canisters/case
10-1090	CaviWipes (220 per can)	12 canisters/case
13-1155	CaviWipesXL (individual) (50 wipes per box)	6 boxes/case
13-1224	CaviWipes in a Flat Pack (45 wipes per pack)	20 packs/case

Composition/Information on Ingredients

Component	CAS No.	Amount
Water	7732-18-5	70-80%
Isopropanol	67-63-0	17.2%
Ethylene Glycol Monobutyl Ether (2-Butoxyethanol)	111-76-2	1-5%
Diisobutylphenoxyethoxyethyldimethylbenzyl ammonium chloride	121-54-0	0.28%

Handling and Storage

Precautions for Safe Handling: Avoid contact with eyes. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse.

Liquid saturant is a combustible liquid. Keep away from heat, open flames and all other sources of ignition. Keep containers closed when not in use. Do not reuse towelette.

Conditions for Safe Storage, Including any Incompatibilities: Store in a cool, well ventilated area away from heat, oxidizers and all sources of ignition.