

CaviWipes1™

1-Minute, 1-Step Disinfecting Towelettes

Fast acting, convenient, ready-to-use, intermediate-level disinfecting wipe.

- Effective against TB, MDR *A. baumannii* and over 20 other microorganisms in just 1 minute*
- 1-Step cleaning and disinfecting allows for faster room turnover
- Low alcohol formula offers excellent materials compatibility
- Stays fully saturated and will not dry out during use
- For use on hard, nonporous surfaces, even in neonatal and infant care areas



Part No.	Product Description	Packaging
13-5100	CaviWipes1 160 count (6" x 6.75")	12 canisters/case
13-5155	CaviWipes1 EXTRA LARGE Singles, 50 count (9" x 12")	6 boxes/case
13-5224	CaviWipes1 Flat Pack, 45 count (7" x 9")	20 packs per case
13-1175	CaviWipes Canister Wall Bracket	12 brackets/case



CaviWipes1 & CaviWipes1
Extra Large can mount on
walls with canister brackets
(13-1175)

Composition/Information On Ingredients

Component	CAS No.	Amount
Water	7732-18-5	70-80%
Isopropanol	67-63-0	15%
Ethanol	64-17-5	7.5%
Ethylene Glycol Monobutyl Ether (2- Butoxyethanol)	111-76-2	1-5%
Didecyldimethylammonium chloride	7173-51-5	0.76%

Handling and Storage

Precautions for Safe Handling: Do not get in eyes or on clothing. Avoid prolonged contact with skin. Use with adequate ventilation. Wash thoroughly with soap and water after handling. Keep containers closed when not in use. Wash contaminated clothing before reuse.

Conditions for Safe Storage, Including any Incompatibilities: Store in a cool, well ventilated area away from heat and all sources of ignition.

* Visit www.metrex.com for a complete list of microbicidal activity