



EmPower™
Dual-Enzymatic Detergent

A dual-enzymatic detergent that works fast to remove blood, tissue, mucous and other protein-rich body fluids from instruments.

- Contains two proteolytic enzymes for effective cleaning
- Gentle on instruments
- Low-foaming formula is ideal for use in automated machines

Part No.	Product Description	Packaging
10-4100	EmPower – 1 gallon	4 bottles/case
10-4150	EmPower – 5 gallon	1 drum

Composition/Information On Ingredients

Component	CAS No.	Amount
Water and non-hazardous ingredients	Mixture	50-60%
Propylene glycol	57-55-6	10-35%
Trisodium Citrate Dihydrate	68-04-2	5-10%
Proteinase subtilisin	9014-01-1	0.1-<1%

Handling and Storage

Precautions for Safe Handling: Avoid contact with eyes, skin and clothing. Avoid generating and breathing mists. Wear appropriate eye protection when handling. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse.

Conditions for Safe Storage, Including any Incompatibilities: Store at room temperature.