



MetriCide™ 28

**High-Level Disinfectant/Sterilant
(2.5% Glutaraldehyde)**

Reusable sterilant and high-level disinfectant for up to 28 days.

- Recommended for immersible and heat-sensitive equipment and instruments
- Contains surfactants

<i>Part No.</i>	<i>Product Description</i>	<i>Packaging</i>
10-2800	MetriCide 28 – 1 gallon	4 bottles/case

Composition/Information on Ingredients

Component	CAS No.	Amount
Water and non-hazardous ingredients	Mixture	97.5%
Glutaraldehyde	111-30-8	2.5%

Handling and Storage

Precautions for Safe Handling: Avoid breathing vapors or mists. Avoid contact with eyes, skin and clothing. Wash thoroughly with soap and water after handling. Avoid contamination of food. Use in a well-ventilated area in closed containers. Use in accordance with label directions.

Conditions for Safe Storage, Including any Incompatibilities: Store at room temperature. Keep containers closed.